


Finally, a urine collection option that can be worn discreetly and emptied while standing!

For Post-Procedure Recovery or Long-Term Catheterization

The Stand Tall™ Urinary Management System gives temporary and long-term Foley catheter users the freedom to empty collected urine while standing at a urinal.

- Thigh location eliminates awkward bending to empty
- Secure waist & leg bands prevent bag migration
- Conceals easily under clothing
- Tapered bag facilitates easy drainage
- Can even be worn unseen with summer shorts

Secure Adjustable Waist Band
Discreet Conceal & Carry™ Leg Bag
Drainage Tube with Activation Valve
Foley Catheter Port with Anti-Reflux Valve
Adjustable Leg Band



Stand and Deliver! The Stand Tall System empties by simply opening the activation valve and putting hand pressure on the bag to drain while **STANDING**.

- No need to remove clothes to empty! -





The Stand Tall Urinary Management System includes the Conceal & Carry™ Leg Bag with all ports and connections ready to use with any standard Foley catheter.

Urologists Prefer “Standing” Drainage

20 out of 20 urologists surveyed preferred a secure leg bag system that allows patients to stand at a urinal or toilet during drainage, compared to existing systems worn on the lower leg.*

* Data on File with United Urologic Innovations, LLC



Visit www.StandTallUrinaryBag.com for pricing and ordering information.
Email: info@StandTallUrinarybag.com



Stand Tall is a trademark of United Urologic Innovations, LLC.
©2018 United Urologic Innovations, LLC
Part No. AA001 Rev B May2018

National Distributor:

NorthStar Surgical
4978 Lincoln Drive
Edina, MN 55436

Toll-Free 866.928.1490
Fax 952.928.1945
www.northstarsurgical.com

